



FCS 500 CEILING CARRIER

HYBRID MONITOR HOLDER













Hillrom offers a variety of arm systems that allow for flexible arrangement and positioning of a various number of monitors. Options range from one large monitor up to eight smaller monitors to fit your OR needs.



	FCS 500-SI	FCS 500-TP	FCS 500-PH
Monitor Compatibility	Siemens	Third Party (excluding Siemens or Philips)	Philips
Support Arm Lengths	1840 mm (72"), 2105 mm (83"), 2370 mm (93"), 2540 mm (100")	1840 mm (72"), 2105 mm (83"), 2370 mm (93"), 2540 mm (100")	1600 mm (63"), 1865 mm (73"), 2130 mm (84"), 2300 mm (91")
Vertically Articulating Arm Type	Spring	Spring	Motorized
Vertically Articulating Arm Length	1,240 mm	1,240 mm	1,000 mm
Maximum Reach	2540 mm (100")	2540 mm (100")	2300 mm (90")
Monitor Configurations	1 - 8 monitors (see back for configuration types)	1 - 8 monitors (see back for configuration types)	1 large monitor (Philips MCS)
Ability to Mount Monitors to Backside	2 small monitors, up to 26"	2 small monitors, up to 26"	2 small monitors, up to 26"
Arm Brake Options <i>*Friction required for monitor bearing brake</i>	Friction or Electro-Pneumatic (gas-assisted spring arm)	Friction or Electro-Pneumatic (gas-assisted spring arm)	Friction, Electro-Pneumatic, (motorized articulating arm)



FCS 500 CONFIGURATIONS

	Type 1	Type 2	Type 3	Type 4	Type 5	Type 6
Siemens Monitor Holders						
Third-Party Monitor Holders						
# of Monitors	3 Monitors	4 Monitors	2 x 2 Monitors	2 x 3 Monitors	2 x 4 Monitors	up to 60" monitor*
Spring Arm Load Capacity	20-50 kg (44-110 lb)	40-70 kg (88-154 lb)	40-70 kg (88-154 lb)	60-90 kg (132-198 lb)	80-110 kg (176-242 lb)	80-110 kg (176-242 lb)
Monitor Carrier Load Capacity	32 kg (max. 70 lb)	36 kg (max. 79 lb)	36 kg (max. 79 lb)	62 kg (max. 136 lb)	62 kg (max. 136 lb)	72 kg (max. 158 lb)

*Must be checked individually

Material Number / REF #	Material Name
#4037261	FCS 500 Ceiling Carrier HL Solo/TanPrep
#4037262	FCS 500 Ceiling Carrier HL TanAdd
#4037251	FCS 500 Ceiling Carrier ML Solo/TanPrep
#4037252	FCS 500 Ceiling Carrier ML TanAdd

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

The FCS Ceiling Supply Unit is intended to be used by clinicians and medically qualified personnel. Class IIb Notified Body: TÜV Rheinland LGA Products GmbH (CE0197), legal manufacturer TRUMPF Medizin Systeme GmbH + Co. KG, Carl-Zeiss-Straße 7-9, 07318 Saalfeld, Germany). This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. TRUMPF Medizin Systeme GmbH + Co. KG is a subsidiary of Hill-Rom Holdings, Inc.