



Hillrom™

Overbed Table - TA270

SAFE, STURDY & ERGONOMIC DESIGN

- The Autolift™ system from Hillrom helps reduce the risk of entrapment and overturning. With a simple push, the table top lifts and moves securely into place.
- The H and U-shaped base deliver maximum stability and are compatible in all care environments.
 - H-shaped design features a reduced width front wheelbase.
 - U-shaped design is ideal for use with a variety of chairs.
- Locking dual-directional tilted table top provides an ideal reading position.
- Plastic tray is designed for easy cleaning and water retention from glass spills.

**High Pressure
Laminate (HPL)
Tray**



Plastic Tray





Hillrom™

TECHNICAL SPECIFICATIONS

Height adjustment	HPL 730 – 1080 mm (without rails) PLASTIC 740 – 1090 mm (H-shaped foot) 752 – 1102 mm (U-shaped foot)
Sizes	Width: 960 mm Depth: 490 mm
Surface area	Fixed models fixed table: 960 x 450 mm Tilting models tilting portion: 740 x 450 mm / fixed portion: 220 x 450 mm
Weight	18.5 kg
Maximum load	16 kg
Material	SURFACE AREA HPL tray – High Pressure Laminate Plastic tray – Polypropylene and UV resistance UNDERCARRIAGE Metallic grey epoxy coated steel (with plastic cover for H-shaped foot)
Casters	4 x 50 mm double-band with 2 brakes
Tilting range	Model TA270 TT max. 30° locked positions / 360° full rotational
Colors	See colors sheet for HPL tray / Grey (RAL 7035) for plastic tray



U-shaped foot



H-shaped foot

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer: Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR86801-EN-EMEA-R1 17-JUN-2020 ENG - EMEA